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## ENGLISH

- Always wear gloves during the insertion procedure. Open the pack and take out the NIO. Make sure that the NIO is free of all packaging parts.
- Select one of the following insertion sites:
  - 2A. Primary site for Intraosseous insertion: Proximal tibia**  
Approximately 1 inch or 2 cm medially and 1/2 inch or 1 cm proximally to the tibial tuberosity.
  - 2B. Secondary site for Intraosseous insertion: Humeral head**  
Adduct the patient’s hand and locate the greater tubercle next to the head of the humerus.  
NOTE: To prevent accidental removal of the device, following the procedure immobilize the patient arm.
- Disinfect the skin on the insertion site by following institutional protocols. Note that this is a two-handed procedure. Place your non-dominant hand on the textured dots located on the lower part of the NIO and position the NIO at a 90 degree angle to the skin at the injection site. The non-dominant hand should maintain this position throughout this procedure.
- Unlock the NIO by rotating the safety cap 90 degree in either direction.
- Place the palm of your dominant hand over the safety cap. Press the device against the patient’s skin and apply downward pressure. While maintaining downward pressure on the device, pull the trigger wings upwards. This action will activate the device.
- Gently pull the NIO up in a rotary motion while holding the base of the needle stabilizer against the insertion site.
- While holding the needle stabilizer base and cannula in place, remove the trocar by pulling it up (use a twisting motion of the trocar if necessary). The trocar removal notch on the distal end of the NIO can be used to assist in removing the trocar from the cannula. Place trocar into an appropriate biohazard container.
- Connect a syringe and confirm secure fitting. If desired, aspirate bone marrow. Always confirm successful needle placement by flushing with up to 20 cc of fluid or per your protocol. It is recommended to use the NIO Fixation sticker to affix the NIO stabilizer. (A). Connect any standard system for infusion (B). Cover the insertion site with a sterile, occlusive dressing.
- Removal instruction: Remove the cannula and needle stabilizer base by twisting and pulling vertically. Dispose of the removed components in the appropriate biohazard container. Cover the insertion site with a sterile, occlusive wound dressing.

### NIO DESIGN

The NIO is an automatic, single use, spring loaded IO device with a double safety mechanism, maximizing caregiver and patient safety. After activation, a unique needle stabilizer fixates the needle firmly.

### INDICATIONS FOR USE

**NIO Adult (NIO-A)** The NIO is intended to provide intraosseous access in the proximal tibia, as an alternative to IV access during emergencies. The device is for use in adult patients only. The NIO is indicated for use in providing intraosseous access as an alternative to IV access during emergencies. Humeral head IO access is indicated when rapid fluid or pharmacological resuscitation is required and intravenous access is not possible. The device is for use in adult patients only.

**RECOMMENDATIONS** • Recommended penetration depth for Prox. Tibia / Humeral head 2.5cm. • Reassess IO site immediately after step 8 (confirm needle placement). It is recommended to frequently monitor the limb every 10 minutes for the first half hour or longer after beginning of drug administration. • Use a pressure bag or blood pressure cuff for optimal infusion rates. • Continue to monitor extremity for complications on a regular basis, especially pre and post infusion. Before drug administration, the IO needle should be checked for placement and patency by confirming the needle is stable in the bone. Bedside Doppler ultrasonography may be useful to ascertain IO placement and flow. • For conscious patients, consider local anesthesia such as lidocaine, per your protocols/policies. • Following IO injection in the Humeral head, in order to prevent accidental removal of the device, immobilize the patient arm. • The intraosseous catheter should be removed within 72 hours (depending on local protocol). Needles should be removed as soon as permanent venous access is established. • Carry backup supplies of critical medical devices whenever possible.

**⚠️ WARNINGS** • Do not aim the NIO toward the joint space or epiphyseal plate. • Discontinue infusion if any signs of infiltration are apparent, including tissue swelling around the insertion site. • The NIO contains sharp parts that should be disposed of in an appropriate container for disposable of medical biohazard waste. • The use of the NIO is restricted to skilled, authorized mediics, nurses, paramedics and doctors who were trained on the device. • Caution: This device is to be used on the order of or by a licensed physician or licensed practitioner. Rx Only. • Do not use if package is damaged. • The device cannot be reused or re-sterilized. Reusing this disposable device might lead to infection, mechanical failure, or harm to the operator and/or patient. • When using any intraosseous device, the possibility of air embolism exists. • Metal needles are not MRI compatibles. • Metal needles may cause scatter artifacts on computed tomography (CT) scans. • The safe use of the NIO in patients with osteoporosis, osteopetrosis, Osgood-Schlatter disease, pre-existing bone pathology, or deformity has not been proven. These conditions may obscure landmarks of the tibia.

**RECOMMENDED STORAGE CONDITIONS** • The NIO device should be stored at room temperature. • The NIO has been tested successfully for storage in an extreme temperature range from (-40) – (+55)°C. • The boxes should be handled with care when placed in a storage room.

**CONTRAINDICATIONS** Stop use of the device should any of the following diagnoses or occurrences arise: • Skin infection at the site location • Tumor • Abnormalities of bone density (e.g. osteogenesis imperfecta, osteopetrosis, osteoporosis) • Osgood-Schlatter disease • Deformation of insertion site • Previous intraosseous insertion / failure on the same bone within last 48h • Previous orthopedic procedures near the insertion site • Fracture of the bone with in the same extremity or the selected bone for insertion • inability to locate anatomical landmarks or excessive tissue

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