



THERMAL PROTECTION? We've Got You Covered!


BLIZZARD - THE ULTIMATE IN HYPOTHERMIA PREVENTION

We've set the standard within the U.S. Military and EMS in thermal performance with a full line of lightweight, reusable weather-proof protection in the most compact of packages.



RADICAL NEW TECHNOLOGY:

Reflexcell™ technology provides maximum warmth by utilizing its cellular construction to trap warm air. Elasticity draws the material to the body and reduces cold spaces and heat loss by convection. Its silvered surfaces block heat loss by radiation. The material is far warmer than goose down, yet completely waterproof and windproof and can be compressed and stored indefinitely without damaging performance.




Our customers include:



ROYAL AIR FORCE

Blizzard also supplies numerous search & rescue teams and Armed forces throughout Europe, Australia & New Zealand.



Applications:

- Medical evacuation
- Casualty treatment
- Personal survival
- Disaster relief
- Search & Rescue



Pre-Hospital Life Support (PHTLS) Military Edition

states that in ALL trauma cases:

"After any immediate life-threatening issues are addressed ... the casualty should be wrapped in a Blizzard Survival Blanket."

The guide also places the Blizzard Survival Blanket FIRST in its 'hierarchical equipment list for prevention and treatment of hypothermia'.

