

Pocket BVM™, PBVM-C

Collapsible Bag Valve Mask • Black Case

- Full 1.6L adult bag valve mask
- Contoured adult mask (size 4) provides a tight seal
- Delivers 500-600cc
- Packed size reduction of 75% compared to standard BVMs
- Can be stored unconventionally
- PEEP compatible



Compact with Impactful Performance

The Pocket BVM™ is a disposable, collapsible, medical-grade silicone resuscitator for adult patients which expands from a disc-like hard shell case to a full-size BVM.

The Pocket BVM reduces cube space by up to 75% over traditional BVMs, is easily stored in unconventional locations and has a ruggedized case that withstands approximately 160lbs psi of direct force. The bag inflates quickly and the mask accommodates most adults. The integrated mask utilizes a curved mold design that mimics the feel and function of inflatable masks without the potential for deflation and delivers a proper seal with each use.



Pocket BVM™, PBVM-C



Product Highlights

Full-size 1.6L bag

Folds compactly into hard shell case


Positive End Expiratory Pressure (PEEP) compatible

This compact, full-size bag valve mask is an ideal component of any medic's kit

- Full size 1.6L adult bag valve mask (standard 15/22 connection)
- Positive End Expiratory Pressure (PEEP) compatible
- Delivers 500-600cc
- Size and weight reduction of 75% compared to traditional BVMs
- Can be stored unconventionally compared to other BVMs
- Contoured silicone adult size 4 mask
- Transparent mask provides a clear view of the procedure
- Hardshell case withstands up to 160lbs of direct force

Specifications

Description	Pocket bag valve mask, black case
Part Number	PBVM-C
NSN	6515-01-590-8909
Packed Size	2.25" x 5.25" (57mm height x 134mm diameter)
Weight	15oz (425g)
Resuscitator Volume	1600ml
Body Mass Range	>40kg (88lbs)
Patient Connector	Inside 15mm female, outside 22mm male
PEEP Connection	30mm male
Reservoir Bag Volume	2600ml
Shelf Life	N/A



PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com
www.ps-med.com

