

T3™ – Next Generation Bandage

Essential and Intuitive Compression

- Multifunctional bandage for various wounds
- Expandable and detachable sterile pad
- Wrap severe wounds and amputations
- Foldable gauze suitable for cleaning and dressing wounds
- Compact, lightweight, waterproof, and vacuum sealed packaging



A Compression Kit in the Palm of Your Hand

The T3™ Bandage offers a compact, multifunctional solution to effectively treat a wide range of injuries including conventional wounds, abdominal and large wounds, amputations, as well as entry and exit wounds.

The T3's versatile design also provides an effective solution for cleaning, covering, and compressing single or multiple wounds.

T3™ – Next Generation Bandage

Product Highlights

Based on the trusted design of the original Emergency Bandage®

Multiple functions in one compact device

Intuitive design facilitates easy, rapid application




This compact, all-in-one bandage is an ideal component of any first aid combat kit.

- Immediate solution for a wide range of injuries
- Easy, intuitive application
- Compact packaging reduces weight and volume
- Integrated pressure applicator applies direct pressure
- Expandable and detachable sterile pad
- Additional gauze for wound cleaning or for use as a secondary dressing
- Elasticized wrapping leader for compression
- Stop and Go Release for controlled wrapping
- Simple closure bar

Specifications

Description	4" (10.16cm) customizable bandage with expandable pad and additional gauze
Part Number	FCP-T3
Pad Size	4" wide x 7" long (10cm x 18cm)
Detachable Pad Size	4" wide x 7" long (10cm x 18cm)
Gauze Size	4" wide x 29.5" long (10cm x 75cm)
Elasticized Leader (expanded)	4" wide x 145" long (10cm x 370cm)
Shelf Life	8 years
Colors	OD Green elastic leader, white pad, white gauze



PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com
www.ps-med.com

